

Smarter Thinking. Simpler Design.



# PrimaConnex<sup>®</sup> & PrimaSolo<sup>®</sup>

Product Catalog





**PrimaConnex®**

Ø3.5 mm PrimaConnex® Implants ..... 6  
 Ø3.5 mm Prosthetic Options ..... 7  
 Ø3.5 mm Prosthetic Components ..... 8  
  
 Ø4.1 mm PrimaConnex® Implants ..... 12  
 Ø4.1 mm Prosthetic Options ..... 13  
 Ø4.1 mm Prosthetic Components ..... 14  
  
 Ø5.0 mm PrimaConnex® Implants ..... 18  
 Ø5.0 mm Prosthetic Options ..... 19  
 Ø5.0 mm Prosthetic Components ..... 20  
  
 LOCATOR® Prosthetic Options ..... 25  
 LOCATOR R-Tx™ Abutment ..... 26  
 LOCATOR® Abutment ..... 27

**PrimaSolo®**

Ø3.0 mm PrimaSolo® Implants ..... 30  
 Ø3.5 mm SD PrimaSolo® Implants ..... 32  
 Ø4.1 mm RD PrimaSolo® Implants ..... 34  
 Ø5.0 mm WD PrimaSolo® Implants ..... 36  
  
 PrimaConnex® Surgical Kit ..... 38  
  
 PrimaConnex® Guided Surgical Kit..... 44  
  
 Surgical Motor.....48  
  
 Terms & Conditions .....50


Product specifications are subject to change without notice.  
 Items illustrated are not to scale.





# PrimaConnex<sup>®</sup>

Smarter Thinking. Simpler Design.




**PrimaConnex® Implants  
Tapered\***

*Catalog #*

*Implant Size*

<b>15413K</b>	Ø3.5 x 10 mm Length
<b>15414K</b>	Ø3.5 x 11.5 mm Length
<b>15415K</b>	Ø3.5 x 13 mm Length
<b>15416K</b>	Ø3.5 x 15 mm Length
<b>15613K</b>	Ø3.5 x 10 mm Length TC
<b>15614K</b>	Ø3.5 x 11.5 mm Length TC
<b>15615K</b>	Ø3.5 x 13 mm Length TC
<b>15616K</b>	Ø3.5 x 15 mm Length TC



**PrimaConnex® Implants  
Straight\***

*Catalog #*

*Implant Size*

<b>15425K</b>	Ø3.3 x 10 mm Length
<b>15426K</b>	Ø3.3 x 11.5 mm Length
<b>15427K</b>	Ø3.3 x 13 mm Length
<b>15428K</b>	Ø3.3 x 15 mm Length
<b>15625K</b>	Ø3.3 x 10 mm Length TC
<b>15626K</b>	Ø3.3 x 11.5 mm Length TC
<b>15627K</b>	Ø3.3 x 13 mm Length TC
<b>15628K</b>	Ø3.3 x 15 mm Length TC



**Cover Screw**

*Catalog #*

*Description*

<b>15300K</b>	Ø3.5 mm
---------------	---------



**Healing Abutments**

*Catalog #*

*Description*


*Flare*

*Cuff Height*

<b>15301K</b>	Ø3.5 mm	Straight	3.0 mm
<b>15302K</b>	Ø3.5 mm	Straight	5.0 mm
<b>15303K</b>	Ø3.5 mm	Straight	7.0 mm
<b>15317K</b>	Ø3.5 mm	4.0 mm	3.0 mm
<b>15318K</b>	Ø3.5 mm	4.0 mm	5.0 mm
<b>15304K</b>	Ø3.5 mm	5.0 mm	3.0 mm
<b>15305K</b>	Ø3.5 mm	5.0 mm	5.0 mm

\*All PrimaConnex® implant sizes are available with either a 1.0 mm machined collar or Total Coverage (TC).

The Ø3.5 mm PrimaConnex® can be restored with straight, Ø4.0 and Ø5.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe abutments.



**Open & Closed Tray Impression Post**



Catalog #	Description	Flare
45130K	Ø3.5 mm	Straight
45136K	Ø3.5 mm	4.0 mm
45131K	Ø3.5 mm	5.0 mm

**Open Tray Impression Post**



Catalog #	Description	Flare
45904K	Ø3.5/3.8 mm	4.0 mm
45906K	Ø3.5/3.8 mm	5.0 mm

**Closed Tray Impression Post**



Catalog #	Description	Flare
45905K	Ø3.5/3.8 mm	4.0 mm
45907K	Ø3.5/3.8 mm	5.0 mm

**Implant Analog**



Catalog #	Description
45140K	Ø3.5 mm, 1-Pack
45143K	Ø3.5 mm, 10-Pack

**Temporary Abutment**




Catalog #	Description	Flare
45118K	Ø3.5 mm, Locking	4.0 mm
45119K	Ø3.5 mm, Non-Locking	4.0 mm
45112K	Ø3.5 mm, Locking	5.0 mm
45113K	Ø3.5 mm, Non-Locking	5.0 mm

**Ti Temporary Abutment**



Catalog #	Description
45075K	Ø3.5 mm, Locking
45076K	Ø3.5 mm, Non-Locking

**Esthetic Contour Ti Abutment**



Catalog #	Description	Flare	Cuff Height
45046K	Ø3.5 mm, Straight	4.0 mm	1.0 mm
45047K	Ø3.5 mm, Straight	4.0 mm	2.0 mm
45048K	Ø3.5 mm, 15° Angled	4.0 mm	1.0 mm
45049K	Ø3.5 mm, 15° Angled	4.0 mm	2.0 mm
45030K	Ø3.5 mm, Straight	5.0 mm	1.0 mm
45031K	Ø3.5 mm, Straight	5.0 mm	3.0 mm
45032K	Ø3.5 mm, 15° Angled	5.0 mm	1.0 mm
45033K	Ø3.5 mm, 15° Angled	5.0 mm	3.0 mm

**UCLA Abutment**



Catalog #	Description	Cuff Height
45100K	Gold/Plastic Sleeve, Locking	1.0 mm
45101K	Gold/Plastic Sleeve, Non-Locking	1.0 mm



	<b>Replacement Final Abutment Screw</b>	<i>Catalog #</i>	<i>Description</i>		
		<b>45060K</b>	1-Pack		
	<b>Laboratory Screw</b>	<i>Catalog #</i>	<i>Description</i>		
		<b>45061K</b>	10-Pack		
	<b>Quick-Packs</b>	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
		<b>45400K</b>	Ø3.5 mm	5.0 mm	1.0 mm
		<b>45401K</b>	Ø3.5 mm	5.0 mm	2.0 mm
		<b>45402K</b>	Ø3.5 mm	5.0 mm	3.0 mm
		<b>45403K</b>	Ø3.5 mm	5.0 mm	4.0 mm
	<b>Quick-Abutment</b>	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
		<b>45420K</b>	Ø3.5 mm	5.0 mm	1.0 mm
		<b>45421K</b>	Ø3.5 mm	5.0 mm	2.0 mm
		<b>45422K</b>	Ø3.5 mm	5.0 mm	3.0 mm
		<b>45423K</b>	Ø3.5 mm	5.0 mm	4.0 mm
	<b>Quick-Abutment Quick-Cap®</b>	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		<b>45432K</b>	Ø3.5 mm, 1-Pack	5.0 mm	
	<b>Quick-Abutment Temporary/Healing Cap</b>	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		<b>45209K</b>	Ø3.5 mm	5.0 mm	
	<b>Quick-Abutment Analog</b>	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		<b>45445K</b>	Ø3.5 mm	5.0 mm	
	<b>Quick-Abutment Plastic Sleeve</b>	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		<b>45215K</b>	Ø3.5 mm, Locking	5.0 mm	
		<b>45218K</b>	Ø3.5 mm, Non-Locking	5.0 mm	
	<b>Quick-Abutment Ceramic Coping</b>	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		<b>45301K</b>	Ø3.5 mm	5.0 mm	



**Reamer Tool for Quick-Abutment Plastic Sleeve**

*Catalog #*

*Description*

**45311K**

Reamer Tool and Centering Pins



**Prepable Abutment**

*Catalog #*

*Description*

**45042K**

Ø3.5 mm




**Ti-Base Abutment (Short)**

*Catalog #*

*Description*

**45806K**

Ø3.5 mm




**Pink Ti-Base Abutment (Short)**

*Catalog #*

*Description*

**45007K**

Ø3.5/3.8 mm




**Ti-Base Abutment (Long)**

*Catalog #*

*Description*

**45901K**

Ø3.5 mm



**Multi Unit Abutment**

*Catalog #*

*Description*

*Cuff Height*

**45312K**

Ø3.5/3.8 mm, Straight

1.0 mm

**45313K**

Ø3.5/3.8 mm, Straight

2.0 mm

**45314K**

Ø3.5/3.8 mm, Straight

3.0 mm

**45321K**

Ø3.5/3.8 mm, Angled, 17°

3.0 mm

**45322K**

Ø3.5/3.8 mm, Angled, 17°

4.0 mm

**45936K**

Ø3.5/3.8 mm, Angled, 30°

3.0 mm

**45937K**

Ø3.5/3.8 mm, Angled, 30°

4.0 mm

**45938K**

Ø3.5/3.8 mm, Angled, 30°

5.0 mm

**Multi Unit Abutment  
Healing Cap**



*Catalog #*

*Description*

**45330K**

Ø3.5/3.8 mm

**Multi Unit Abutment  
Impression Post**



*Catalog #*

*Description*

**45327K**  
**45328K**

Open Tray, Ø3.5/3.8 mm with screw  
Closed Tray, Ø3.5/3.8 mm

**Multi Unit Abutment  
Analog**



*Catalog #*

*Description*

**45329K**

Ø3.5/3.8 mm

**Multi-Unit Abutment  
Ti Temporary Sleeve**



*Catalog #*

*Description*

**45331K**

Ø3.5/3.8 mm

**Multi Unit Abutment  
Plastic Sleeve**



*Catalog #*

*Description*

**45332K**

Ø3.5/3.8 mm

**Multi Unit Abutment  
Prosthetic Screw**



*Catalog #*

*Description*

**45333K**

1-Pack

**Replacement Abutment  
Screw for Multi Unit  
Angled Abutment**




*Catalog #*

*Description*

**45334K**

1-Pack



**PrimaConnex® Implants  
Tapered\***

*Catalog #*


**15417K**  
**15418K**  
**15419K**  
**15420K**

*Implant Size*

Ø4.1 x 10 mm Length  
Ø4.1 x 11.5 mm Length  
Ø4.1 x 13 mm Length  
Ø4.1 x 15 mm Length

**15617K**  
**15618K**  
**15619K**  
**15620K**

Ø4.1 x 10 mm Length TC  
Ø4.1 x 11.5 mm Length TC  
Ø4.1 x 13 mm Length TC  
Ø4.1 x 15 mm Length TC



**PrimaConnex® Implants  
Straight\***

*Catalog #*

**15429K**  
**15430K**  
**15431K**  
**15432K**  
**15433K**

*Implant Size*

Ø4.0 x 8 mm Length  
Ø4.0 x 10 mm Length  
Ø4.0 x 11.5 mm Length  
Ø4.0 x 13 mm Length  
Ø4.0 x 15 mm Length

**15629K**  
**15630K**  
**15631K**  
**15632K**  
**15633K**

Ø4.0 x 8 mm Length TC  
Ø4.0 x 10 mm Length TC  
Ø4.0 x 11.5 mm Length TC  
Ø4.0 x 13 mm Length TC  
Ø4.0 x 15 mm Length TC



**Cover Screw**

*Catalog #*

**15306K**

*Description*

Ø4.1 mm



**Healing Abutments**

*Catalog #*

**15307K**  
**15308K**  
**15309K**

*Description*

Ø4.1 mm  
Ø4.1 mm  
Ø4.1 mm

*Flare*

Straight  
Straight  
Straight

*Cuff Height*

3.0 mm  
5.0 mm  
7.0 mm

**15319K**  
**15320K**

Ø4.1 mm  
Ø4.1 mm

5.0 mm  
5.0 mm

3.0 mm  
5.0 mm

**15310K**  
**15311K**


Ø4.1 mm  
Ø4.1 mm

6.0 mm  
6.0 mm

3.0 mm  
5.0 mm

\*All PrimaConnex® implant sizes are available with either a 1.0 mm machined collar or Total Coverage (TC).

The Ø4.1 mm PrimaConnex® can be restored with straight, Ø5.0 and Ø6.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe abutments.



**Open & Closed Tray  
Impression Post**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45132K</b>	Ø4.1 mm	Straight
<b>45137K</b>	Ø4.1 mm	5.0 mm
<b>45133K</b>	Ø4.1 mm	6.0 mm

**Open Tray  
Impression Post**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45908K</b>	Ø4.1/4.5 mm	5.0 mm
<b>45910K</b>	Ø4.1/4.5 mm	6.0 mm

**Closed Tray  
Impression Post**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45909K</b>	Ø4.1/4.5 mm	5.0 mm
<b>45911K</b>	Ø4.1/4.5 mm	6.0 mm

**Implant Analog**



<i>Catalog #</i>	<i>Description</i>
<b>45141K</b>	Ø4.1 mm, 1-Pack
<b>45144K</b>	Ø4.1 mm, 10-Pack

**Temporary Abutment**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45120K</b>	Ø4.1 mm, Locking	5.0 mm
<b>45121K</b>	Ø4.1 mm, Non-Locking	5.0 mm
<b>45114K</b>	Ø4.1 mm, Locking	6.0 mm
<b>45115K</b>	Ø4.1 mm, Non-Locking	6.0 mm

**Ti Temporary Abutment**



<i>Catalog #</i>	<i>Description</i>
<b>45077K</b>	Ø4.1 mm, Locking
<b>45078K</b>	Ø4.1 mm, Non-Locking

**Esthetic Contour Ti  
Abutment**



Cuff Height  
(Buccal)

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
<b>45050K</b>	Ø4.1 mm, Straight	5.0 mm	0.5 mm
<b>45051K</b>	Ø4.1 mm, Straight	5.0 mm	1.5 mm
<b>45052K</b>	Ø4.1 mm, 15° Angled	5.0 mm	0.5 mm
<b>45053K</b>	Ø4.1 mm, 15° Angled	5.0 mm	1.5 mm
<b>45034K</b>	Ø4.1 mm, Straight	6.0 mm	1.0 mm
<b>45035K</b>	Ø4.1 mm, Straight	6.0 mm	3.0 mm
<b>45036K</b>	Ø4.1 mm, 15° Angled	6.0 mm	1.0 mm
<b>45037K</b>	Ø4.1 mm, 15° Angled	6.0 mm	3.0 mm

**UCLA Abutment**



<i>Catalog #</i>	<i>Description</i>	<i>Cuff Height</i>
<b>45102K</b>	Gold/Plastic Sleeve, Locking	1.0 mm
<b>45103K</b>	Gold/Plastic Sleeve, Non-Locking	1.0 mm



**Replacement  
Final Abutment Screw**

<i>Catalog #</i>	<i>Description</i>
<b>45060K</b>	1-Pack




**Laboratory Screw**

<i>Catalog #</i>	<i>Description</i>
<b>45061K</b>	10-Pack



**Quick-Packs**

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
<b>45412K</b>	Ø4.1 mm	5.0 mm	0.5 mm
<b>45413K</b>	Ø4.1 mm	5.0 mm	1.5 mm
<b>45414K</b>	Ø4.1 mm	5.0 mm	2.5 mm
<b>45415K</b>	Ø4.1 mm	5.0 mm	4.0 mm
<b>45404K</b>	Ø4.1 mm	6.0 mm	0.5 mm
<b>45405K</b>	Ø4.1 mm	6.0 mm	1.5 mm
<b>45406K</b>	Ø4.1 mm	6.0 mm	2.5 mm
<b>45407K</b>	Ø4.1 mm	6.0 mm	4.0 mm



**Quick-Abutment**

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
<b>45435K</b>	Ø4.1 mm	5.0 mm	0.5 mm
<b>45436K</b>	Ø4.1 mm	5.0 mm	1.5 mm
<b>45437K</b>	Ø4.1 mm	5.0 mm	2.5 mm
<b>45438K</b>	Ø4.1 mm	5.0 mm	4.0 mm
<b>45424K</b>	Ø4.1 mm	6.0 mm	0.5 mm
<b>45425K</b>	Ø4.1 mm	6.0 mm	1.5 mm
<b>45426K</b>	Ø4.1 mm	6.0 mm	2.5 mm
<b>45427K</b>	Ø4.1 mm	6.0 mm	4.0 mm



**Quick-Abutment  
Quick-Cap®**

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45443K</b>	Ø4.1 mm, 1-Pack	5.0 mm
<b>45433K</b>	Ø4.1 mm, 1-Pack	6.0 mm



**Quick-Abutment  
Temporary/Healing Cap**

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45236K</b>	Ø4.1 mm	5.0 mm
<b>45210K</b>	Ø4.1 mm	6.0 mm



**Quick-Abutment  
Analog**

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45448K</b>	Ø4.1 mm	5.0 mm
<b>45446K</b>	Ø4.1 mm	6.0 mm



**Quick-Abutment Plastic Sleeve**

*Catalog #*

*Description*

*Flare*

**45240K**  
**45242K**

Ø4.1 mm, Locking  
Ø4.1 mm, Non-Locking

5.0 mm  
5.0 mm

**45216K**  
**45219K**

Ø4.1 mm, Locking  
Ø4.1 mm, Non-Locking

6.0 mm  
6.0 mm



**Quick-Abutment Ceramic Coping**

*Catalog #*

*Description*

*Flare*

**45244K**  
**45302K**

Ø4.1 mm  
Ø4.1 mm

5.0 mm  
6.0 mm



**Reamer Tool for Quick-Abutment Plastic Sleeve**

*Catalog #*

*Description*

**45311K**

Reamer Tool and Centering Pins



**Prepable Abutment**

*Catalog #*

*Description*

**45043K**

Ø4.1 mm




**Ti-Base Abutment (Short)**

*Catalog #*

*Description*

**45807K**

Ø4.1 mm



**Pink Ti-Base Abutment (Short)**

*Catalog #*

*Description*

**45008K**

Ø4.1/4.5 mm



**Ti-Base Abutment (Long)**

*Catalog #*


*Description*

**45902K**

Ø4.1 mm



**Multi Unit Abutment**



<i>Catalog #</i>	<i>Description</i>	<i>Cuff Height</i>
<b>45315K</b>	Ø4.1/4.5 mm, Straight	1.0 mm
<b>45316K</b>	Ø4.1/4.5 mm, Straight	2.0 mm
<b>45317K</b>	Ø4.1/4.5 mm, Straight	3.0 mm
<b>45318K</b>	Ø4.1/4.5 mm, Straight	4.0 mm
<b>45323K</b>	Ø4.1/4.5 mm, Angled, 17°	3.0 mm
<b>45324K</b>	Ø4.1/4.5 mm, Angled, 17°	4.0 mm
<b>45325K</b>	Ø4.1/4.5 mm, Angled, 30°	4.0 mm
<b>45326K</b>	Ø4.1/4.5 mm, Angled, 30°	5.0 mm

**Multi Unit Abutment Healing Cap**



<i>Catalog #</i>	<i>Description</i>
<b>45330K</b>	Ø4.1/4.5 mm

**Multi Unit Abutment Impression Post**



<i>Catalog #</i>	<i>Description</i>
<b>45327K</b>	Open Tray, Ø4.1/4.5 mm with screw
<b>45328K</b>	Closed Tray, Ø4.1/4.5 mm

**Multi Unit Abutment Analog**



<i>Catalog #</i>	<i>Description</i>
<b>45329K</b>	Ø4.1/4.5 mm

**Multi-Unit Abutment Ti Temporary Sleeve**



<i>Catalog #</i>	<i>Description</i>
<b>45331K</b>	Ø4.1/4.5 mm

**Multi Unit Abutment Plastic Sleeve**



<i>Catalog #</i>	<i>Description</i>
<b>45332K</b>	Ø4.1/4.5 mm

**Multi Unit Abutment Prosthetic Screw**




<i>Catalog #</i>	<i>Description</i>
<b>45333K</b>	1-Pack

**Replacement Abutment Screw for Multi Unit Angled Abutment**



<i>Catalog #</i>	<i>Description</i>
<b>45334K</b>	1-Pack



**PrimaConnex® Implants  
Tapered\***

*Catalog #*


*Implant Size*

<b>15421K</b>	Ø5.0 x 10 mm Length
<b>15422K</b>	Ø5.0 x 11.5 mm Length
<b>15423K</b>	Ø5.0 x 13 mm Length
<b>15424K</b>	Ø5.0 x 15 mm Length
<b>15621K</b>	Ø5.0 x 10 mm Length TC
<b>15622K</b>	Ø5.0 x 11.5 mm Length TC
<b>15623K</b>	Ø5.0 x 13 mm Length TC
<b>15624K</b>	Ø5.0 x 15 mm Length TC

**PrimaConnex® Implants  
Straight\***

*Catalog #*

*Implant Size*



<b>15434K</b>	Ø5.0 x 8 mm Length
<b>15435K</b>	Ø5.0 x 10 mm Length
<b>15436K</b>	Ø5.0 x 11.5 mm Length
<b>15437K</b>	Ø5.0 x 13 mm Length
<b>15634K</b>	Ø5.0 x 8 mm Length TC
<b>15635K</b>	Ø5.0 x 10 mm Length TC
<b>15636K</b>	Ø5.0 x 11.5 mm Length TC
<b>15637K</b>	Ø5.0 x 13 mm Length TC

**Cover Screw**

*Catalog #*

*Description*



<b>15312K</b>	Ø5.0 mm
---------------	---------

**Healing Abutments**

*Catalog #*

*Description*

*Flare*


*Cuff Height*



<b>15313K</b>	Ø5.0 mm	Straight	3.0 mm
<b>15314K</b>	Ø5.0 mm	Straight	5.0 mm
<b>15321K</b>	Ø5.0 mm	5.0 mm	3.0 mm
<b>15322K</b>	Ø5.0 mm	5.0 mm	5.0 mm
<b>15315K</b>	Ø5.0 mm	6.0 mm	3.0 mm
<b>15316K</b>	Ø5.0 mm	6.0 mm	5.0 mm

\*All PrimaConnex® implant sizes are available with either a 1.0 mm machined collar or Total Coverage (TC).

The Ø5.0 mm PrimaConnex® can be restored with straight, Ø5.0 and Ø6.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe abutments.





**Open & Closed Tray  
Impression Post**

*Catalog #*

*Description*

*Flare*

**45134K**  
**45138K**  
**45135K**

Ø5.0 mm  
Ø5.0 mm  
Ø5.0 mm

Straight  
5.0 mm  
6.0 mm



**Open Tray  
Impression Post**

*Catalog #*

*Description*

*Flare*

**45916K**  
**45918K**

Ø5.5/6.5 mm  
Ø5.5/6.5 mm

7.0 mm  
8.0 mm



**Closed Tray  
Impression Post**

*Catalog #*

*Description*

*Flare*

**45917K**  
**45919K**

Ø5.5/6.5 mm  
Ø5.5/6.5 mm

7.0 mm  
8.0 mm



**Implant Analog**

*Catalog #*

*Description*

**45142K**  
**45145K**

Ø5.0 mm, 1-Pack  
Ø5.0 mm, 10-Pack



**Temporary Abutment**

*Catalog #*

*Description*

*Flare*

**45122K**  
**45123K**

Ø5.0 mm, Locking  
Ø5.0 mm, Non-Locking

5.0 mm  
5.0 mm

**45116K**  
**45117K**

Ø5.0 mm, Locking  
Ø5.0 mm, Non-Locking

6.0 mm  
6.0 mm



**Ti Temporary Abutment**

*Catalog #*

*Description*

**45079K**  
**45080K**

Ø5.0 mm, Locking  
Ø5.0 mm, Non-Locking

**Esthetic Contour Ti  
Abutment**

*Catalog #*

*Description*

*Flare*

*Cuff Height*

Cuff Height  
(Buccal)



**45054K**  
**45055K**  
**45056K**  
**45057K**

Ø5.0 mm, Straight  
Ø5.0 mm, Straight  
Ø5.0 mm, 15° Angled  
Ø5.0 mm, 15° Angled

5.0 mm 0.5 mm  
5.0 mm 1.5 mm  
5.0 mm 0.5 mm  
5.0 mm 1.5 mm

**45038K**  
**45039K**  
**45040K**

Ø5.0 mm, Straight  
Ø5.0 mm, Straight  
Ø5.0 mm, 15° Angled

6.0 mm 1.0 mm  
6.0 mm 3.0 mm  
6.0 mm 1.0 mm



**UCLA Abutment**

*Catalog #*

*Description*

*Cuff Height*

**45104K**  
**45105K**

Gold/Plastic Sleeve, Locking  
Gold/Plastic Sleeve, Non-Locking

1.0 mm  
1.0 mm



**Replacement  
Final Abutment Screw**

<i>Catalog #</i>	<i>Description</i>
<b>45060K</b>	1-Pack



**Laboratory Screw**


<i>Catalog #</i>	<i>Description</i>
<b>45061K</b>	10-Pack

**Quick-Packs**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
<b>45416K</b>	Ø5.0 mm	5.0 mm	0.5 mm
<b>45417K</b>	Ø5.0 mm	5.0 mm	1.5 mm
<b>45418K</b>	Ø5.0 mm	5.0 mm	2.5 mm
<b>45419K</b>	Ø5.0 mm	5.0 mm	4.0 mm
<b>45408K</b>	Ø5.0 mm	6.0 mm	0.5 mm
<b>45409K</b>	Ø5.0 mm	6.0 mm	1.5 mm
<b>45410K</b>	Ø5.0 mm	6.0 mm	2.5 mm
<b>45411K</b>	Ø5.0 mm	6.0 mm	4.0 mm

**Quick-Abutment**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
<b>45439K</b>	Ø5.0 mm	5.0 mm	0.5 mm
<b>45440K</b>	Ø5.0 mm	5.0 mm	1.5 mm
<b>45441K</b>	Ø5.0 mm	5.0 mm	2.5 mm
<b>45442K</b>	Ø5.0 mm	5.0 mm	4.0 mm
<b>45428K</b>	Ø5.0 mm	6.0 mm	0.5 mm
<b>45429K</b>	Ø5.0 mm	6.0 mm	1.5 mm
<b>45430K</b>	Ø5.0 mm	6.0 mm	2.5 mm
<b>45431K</b>	Ø5.0 mm	6.0 mm	4.0 mm

**Quick-Abutment  
Quick-Cap®**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45444K</b>	Ø5.0 mm, 1-Pack	5.0 mm
<b>45434K</b>	Ø5.0 mm, 1-Pack	6.0 mm

**Quick-Abutment  
Temporary/Healing Cap**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45237K</b>	Ø5.0 mm	5.0 mm
<b>45211K</b>	Ø5.0 mm	6.0 mm

**Quick-Abutment  
Analog**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
<b>45449K</b>	Ø5.0 mm	5.0 mm
<b>45447K</b>	Ø5.0 mm	6.0 mm



**Quick-Abutment Plastic Sleeve**

*Catalog #*

*Description*

*Flare*

**45241K**  
**45243K**  
**45217K**  
**45220K**

Ø5.0 mm, Locking  
Ø5.0 mm, Non-Locking  
Ø5.0 mm, Locking  
Ø5.0 mm, Non-Locking

5.0 mm  
5.0 mm  
6.0 mm  
6.0 mm



**Quick-Abutment Ceramic Coping**

*Catalog #*

*Description*

*Flare*

**45245K**  
**45303K**

Ø5.0 mm  
Ø5.0 mm

5.0 mm  
6.0 mm



**Reamer Tool for Quick-Abutment Plastic Sleeve**

*Catalog #*

*Description*

**45311K**

Reamer Tool and Centering Pins




**Prepable Abutment**

*Catalog #*

*Description*

**45044K**

Ø5.0 mm




**Ti-Base Abutment (Short)**

*Catalog #*

*Description*

**45808K**

Ø5.0 mm




**Pink Ti-Base Abutment (Short)**

*Catalog #*

*Description*

**45009K**

Ø5.0/5.5/6.5 mm



**Ti-Base Abutment (Long)**


*Catalog #*

*Description*

**45903K**

Ø5.0 mm

**Multi Unit Abutment**



<i>Catalog #</i>	<i>Description</i>	<i>Cuff Height</i>
<b>45319K</b>	Ø5.0/5.5/6.5 mm, Straight	1.0 mm
<b>45320K</b>	Ø5.0/5.5/6.5 mm, Straight	2.0 mm
<b>45939K</b>	Ø5.0/5.5/6.5 mm, Straight	3.0 mm
<b>45940K</b>	Ø5.0/5.5/6.5 mm, Straight	4.0 mm
<b>45941K*</b>	Ø5.0/5.5/6.5 mm, Angled, 17°	3.0 mm
<b>45942K*</b>	Ø5.0/5.5/6.5 mm, Angled, 17°	4.0 mm

**Multi Unit Abutment Healing Cap**



<i>Catalog #</i>	<i>Description</i>
<b>45330K</b>	Ø5.0/5.5/6.5 mm

**Multi Unit Abutment Impression Post**



<i>Catalog #</i>	<i>Description</i>
<b>45327K</b>	Open Tray, Ø5.0/5.5/6.5 mm with screw
<b>45328K</b>	Closed Tray, Ø5.0/5.5/6.5 mm

**Multi Unit Abutment Analog**



<i>Catalog #</i>	<i>Description</i>
<b>45329K</b>	Ø5.0/5.5/6.5 mm

**Multi-Unit Abutment Ti Temporary Sleeve**



<i>Catalog #</i>	<i>Description</i>
<b>45331K</b>	Ø5.0/5.5/6.5 mm

**Multi Unit Abutment Plastic Sleeve**



<i>Catalog #</i>	<i>Description</i>
<b>45332K</b>	Ø5.0/5.5/6.5 mm

**Multi Unit Abutment Prosthetic Screw**



<i>Catalog #</i>	<i>Description</i>
<b>45333K</b>	1-Pack

**Replacement Abutment Screw for Multi Unit Angled Abutment**




<i>Catalog #</i>	<i>Description</i>
<b>45334K</b>	1-Pack

\*Available mid-2018







3.5

LOCATOR® Abutment

4.1

5.0

Internal TiLobe



### LOCATOR R-Tx™ Abutment

#### Catalog #

**LOCR-TX 31200-01**  
**LOCR-TX 31200-02**  
**LOCR-TX 31200-03**  
**LOCR-TX 31200-04**  
**LOCR-TX 31200-05**  
**LOCR-TX 31200-06**

#### Description

Ø3.5/3.8 mm  
 Ø3.5/3.8 mm  
 Ø3.5/3.8 mm  
 Ø3.5/3.8 mm  
 Ø3.5/3.8 mm  
 Ø3.5/3.8 mm

#### Cuff Height

1.0 mm  
 2.0 mm  
 3.0 mm  
 4.0 mm  
 5.0 mm  
 6.0 mm

**LOCR-TX 31201-01**  
**LOCR-TX 31201-02**  
**LOCR-TX 31201-03**  
**LOCR-TX 31201-04**  
**LOCR-TX 31201-05**  
**LOCR-TX 31201-06**

Ø4.1/4.5 mm  
 Ø4.1/4.5 mm  
 Ø4.1/4.5 mm  
 Ø4.1/4.5 mm  
 Ø4.1/4.5 mm  
 Ø4.1/4.5 mm

1.0 mm  
 2.0 mm  
 3.0 mm  
 4.0 mm  
 5.0 mm  
 6.0 mm

**LOCR-TX 31202-01**  
**LOCR-TX 31202-02**  
**LOCR-TX 31202-03**  
**LOCR-TX 31202-04**  
**LOCR-TX 31202-05**  
**LOCR-TX 31202-06**

Ø5.0/5.5/6.5 mm  
 Ø5.0/5.5/6.5 mm  
 Ø5.0/5.5/6.5 mm  
 Ø5.0/5.5/6.5 mm  
 Ø5.0/5.5/6.5 mm  
 Ø5.0/5.5/6.5 mm

1.0 mm  
 2.0 mm  
 3.0 mm  
 4.0 mm  
 5.0 mm  
 6.0 mm

### LOCATOR R-Tx™ Abutment Accessories



#### Catalog #

**LOCR-TX 30013-01**

Denture Attachment Processing Assembly, 4-Pack

**LOCR-TX 30001-01**

Zero Retention Insert, 4-Pack

**LOCR-TX 30002-01**

Low Retention Insert, 4-Pack

**LOCR-TX 30003-01**

Medium Retention Insert, 4-Pack

**LOCR-TX 30004-01**

High Retention Insert, 4-Pack

**LOCR-TX 30005-01**

Extra-High Retention Insert, 4-Pack

**LOCR-TX 30014-01**

3.35 mm Abutment Analog, 4-Pack

**LOCR-TX 30015-01**

4.0 mm Abutment Analog, 4-Pack

**LOCR-TX 30016-01**

5.0 mm Abutment Analog, 4-Pack

**LOCR-TX 30017-01**

Impression Coping, 4-Pack

**LOCR-TX 30021-01**

Retention Insert Tool

**LOCR-TX 04936**

.050"/1.25 mm Hex, Torque Wrench

**LOCR-TX 08007**

.050"/1.25 mm Hex Driver with Thumb Knob

**LOCR-TX 08008**

.050"/1.25 mm Hex, Latch Lock Driver

**LOCR-TX 30012-01**

Processing Insert, 4-Pack

**LOCR-TX 30018-01**

Processing Spacer, 4-Pack



**LOCATOR® Abutment Assembly**

<i>Catalog #</i>	<i>Description</i>	<i>Cuff Height</i>
<b>45062K</b>	Ø3.5 mm	1.0 mm
<b>45258K</b>	Ø3.5 mm	2.0 mm
<b>45259K</b>	Ø3.5 mm	3.0 mm
<b>45260K</b>	Ø3.5 mm	4.0 mm
<b>45063K</b>	Ø3.5 mm	5.0 mm
<b>45064K</b>	Ø3.5 mm	6.0 mm
<b>45065K</b>	Ø4.1 mm	1.0 mm
<b>45261K</b>	Ø4.1 mm	2.0 mm
<b>45262K</b>	Ø4.1 mm	3.0 mm
<b>45263K</b>	Ø4.1 mm	4.0 mm
<b>45066K</b>	Ø4.1 mm	5.0 mm
<b>45067K</b>	Ø4.1 mm	6.0 mm
<b>45068K</b>	Ø5.0 mm	1.0 mm
<b>45264K</b>	Ø5.0 mm	2.0 mm
<b>45265K</b>	Ø5.0 mm	3.0 mm
<b>45069K</b>	Ø5.0 mm	4.0 mm
<b>45070K</b>	Ø5.0 mm	5.0 mm
<b>45071K</b>	Ø5.0 mm	6.0 mm

**LOCATOR® Abutment Accessories**




<i>Catalog #</i>	<i>Description</i>
<b>L85192</b>	Denture Cap with Black Processing Male, 2-Pack
<b>45266K</b>	LOCATOR® Replacement Male, White, Standard Retention, 4-Pack
<b>45267K</b>	LOCATOR® Replacement Male, Pink, Light Retention, 4-Pack
<b>45268K</b>	LOCATOR® Replacement Male, Blue, (ELR) Extra Light Retention, 4-Pack
<b>45269K</b>	LOCATOR® Replacement Male, Green, (ER) Extended Range, 4-Pack
<b>45270K</b>	LOCATOR® Female Analog, 4-Pack
<b>45271K</b>	LOCATOR® Impression Coping, 4-Pack
<b>45273K</b>	LOCATOR® Core Tool
<b>45274K</b>	LOCATOR® Torque Wrench Driver





# PrimaSolo®



 <p>Ø3.0 mm</p> <p>2.3 mm</p>	<b>PrimaSolo® Implants Tapered*</b>	<i>Catalog #</i>	<i>Implant Size</i>
		<b>15400K</b>	Ø3.0 x 13 mm Length
		<b>15401K</b>	Ø3.0 x 15 mm Length



<b>Temporary/Healing Cap</b>	<i>Catalog #</i>	<i>Description</i>
	<b>45000K</b>	Ø3.0 mm



<b>Implant Analog</b>	<i>Catalog #</i>	<i>Description</i>
	<b>45172K</b>	Ø3.0 mm

\*PrimaSolo® One-Piece Implants consist of a threaded, tapered implant and integrated gold-colored abutment, placed in a single-stage procedure.

Angle correction up to 10° for Ø3.0 PrimaSolo® Implants; up to 15° for Ø3.5, Ø4.1, and Ø5.0 PrimaSolo® Implants.

Implant Analog



Temporary/Healing Cap



Ø3.0 Implant





Ø3.5 mm



2.4 mm

**PrimaSolo® Implants  
Tapered\***

*Catalog #*

**15402K**

**15403K**

**15404K**

**15405K**

*Implant Size*

Ø3.5 x 10 mm Length

Ø3.5 x 11.5 mm Length

Ø3.5 x 13 mm Length

Ø3.5 x 15 mm Length



**Temporary/Healing Cap**

*Catalog #*

**45001K**

*Description*

Ø3.5 mm



**Implant Analog**

*Catalog #*

**45173K**

*Description*

Ø3.5 mm

\*PrimaSolo® One-Piece Implants consist of a threaded, tapered implant and integrated gold-colored abutment, placed in a single-stage procedure.

Angle correction up to 15° for Ø3.5, Ø4.1, and Ø5.0 PrimaSolo® Implants.



Implant Analog



Temporary/Healing Cap



Ø3.5 Implant





**PrimaSolo® Implants  
Tapered\***

*Catalog #*

*Implant Size*

**15406K**

Ø4.1 x 10 mm Length

**15407K**

Ø4.1 x 11.5 mm Length

**15408K**

Ø4.1 x 13 mm Length

**15409K**

Ø4.1 x 15 mm Length



**Temporary/Healing Cap**

*Catalog #*

*Description*

**45002K**

Ø4.1 mm



**Implant Analog**

*Catalog #*

*Description*

**45174K**

Ø4.1 mm

\*PrimaSolo® One-Piece Implants consist of a threaded, tapered implant and integrated gold-colored abutment, placed in a single-stage procedure.

Angle correction up to 15° for Ø3.5, Ø4.1, and Ø5.0 PrimaSolo® Implants.

Implant Analog



Temporary/Healing Cap



Ø4.1 Implant





**PrimaSolo® Implants  
Tapered\***

*Catalog #*

*Implant Size*

**15410K**

Ø5.0 x 10 mm Length

**15411K**

Ø5.0 x 11.5 mm Length

**15412K**

Ø5.0 x 13 mm Length



**Temporary/Healing Cap**

*Catalog #*

*Description*

**45003K**

Ø5.0 mm



**Implant Analog**

*Catalog #*

*Description*

**45175K**

Ø5.0 mm

\*PrimaSolo® One-Piece Implants consist of a threaded, tapered implant and integrated gold-colored abutment, placed in a single-stage procedure.

Angle correction up to 15° for Ø3.5, Ø4.1, and Ø5.0 PrimaSolo® Implants.

Implant Analog




Temporary/Healing Cap



Ø5.0 Implant





*Catalog # Surgical Options*

<b>15278K</b>	Surgical Cassette (empty)
<b>15241K</b>	<b>A.</b> PrimaConnex® Surgical Kit
<b>15244K</b>	<b>B.</b> PrimaConnex® Tapered Drill and Tap Set
<b>15249K</b>	<b>C.</b> PrimaConnex® Straight Drill and Tap Set

The PrimaConnex® Surgical Kit contains all instrumentation except the final tapered and/or straight drills and taps. The Tapered and/or Straight Drill and Tap Set(s) must be ordered separately, depending on the implant design to be placed.

*Catalog # Description*

<b>15241K</b>	<b>A.</b> PrimaConnex® Surgical Kit <i>Includes:</i>
<b>15278K</b>	Surgical Cassette
<b>15200K</b>	1.5 x 15 mm Initial Drill
<b>15201K</b>	1.8 x 11.5 mm Initial Drill
<b>15202K</b>	1.8 x 15 mm Initial Drill
<b>15225K</b>	Depth Probe
<b>L1404-01K</b>	Drill Extender
<b>15286K</b>	3.5/3.8 Implant Driver, TiLobe Connection, Short
<b>15287K</b>	3.5/3.8 Implant Driver, TiLobe Connection, Long
<b>15288K</b>	4.1/4.5 Implant Driver, TiLobe Connection, Short
<b>15289K</b>	4.1/4.5 Implant Driver, TiLobe Connection, Long
<b>15290K</b>	5.0/5.5/6.5/7.0/8.0/9.0 Implant Driver, TiLobe Connection, Short
<b>15282K</b>	3.5 Implant Driver, PrimaConnex®, Ratchet
<b>15283K</b>	4.1 Implant Driver, PrimaConnex®, Ratchet
<b>15284K</b>	5.0 Implant Driver, PrimaConnex®, Ratchet

*Catalog # Description*


<b>L1800-22-02K</b>	Irrigation Port Cleaner
<b>15226K</b>	Parallel Pins, Straight, 1.5 mm/1.8 mm, 2-Pack
<b>45045K</b>	Parallel Pins, 15° Angled, 1.5 mm/1.8 mm, 2-Pack
<b>15274K</b>	Quad Driver, Latch, Short
<b>15275K</b>	Quad Driver, Latch, Long
<b>15272K</b>	Quad Driver, Swivel, Short
<b>15273K</b>	Quad Driver, Swivel, Long
<b>15224K</b>	Ratchet Adapter
<b>L1702-01K</b>	Surgical Ratchet
<b>L1407K</b>	Conversion Handle, Latch
<b>4409-0255</b>	PrimaConnex® Surgical Manual
<b>15240K</b>	Implant Sizing Overlays
<b>4409-0257</b>	Wall Chart
<b>4409-0258</b>	Drill Usage Chart



Catalog #	Description
<b>15244K</b>	<b>B. PrimaConnex® Tapered Drill and Tap Set</b> <i>Includes:</i>
<b>15205K</b>	3.5 x 10 mm Final Tapered Drill
<b>15206K</b>	3.5 x 11.5 mm Final Tapered Drill
<b>15207K</b>	3.5 x 13 mm Final Tapered Drill
<b>15208K</b>	3.5 x 15 mm Final Tapered Drill
<b>15209K</b>	4.1 x 10 mm Final Tapered Drill
<b>15210K</b>	4.1 x 11.5 mm Final Tapered Drill
<b>15211K</b>	4.1 x 13 mm Final Tapered Drill
<b>15212K</b>	4.1 x 15 mm Final Tapered Drill
<b>15213K</b>	5.0 x 10 mm Final Tapered Drill
<b>15214K</b>	5.0 x 11.5 mm Final Tapered Drill
<b>15215K</b>	5.0 x 13 mm Final Tapered Drill
<b>15075K</b>	5.0 x 15 mm Final Tapered Drill
<b>15217K</b>	3.5 x 10/11.5 mm Tapered Surgical Tap
<b>15247K</b>	3.5 x 13/15 mm Tapered Surgical Tap
<b>15218K</b>	4.1 x 10/11.5 mm Tapered Surgical Tap
<b>15248K</b>	4.1 x 13/15 mm Tapered Surgical Tap
<b>15219K</b>	5.0 x 10/11.5 mm Tapered Surgical Tap
<b>15076K</b>	5.0 x 13/15 mm Tapered Surgical Tap

Catalog #	Description
<b>15249K</b>	<b>C. PrimaConnex® Straight Drill and Tap Set</b> <i>Includes:</i>
<b>15250K</b>	2.8 mm x 11.5 mm Straight Drill
<b>15251K</b>	2.8 mm x 15 mm Straight Drill
<b>15252K</b>	3.3 mm x 11.5 mm Straight Drill
<b>15253K</b>	3.3 mm x 15 mm Straight Drill
<b>15254K</b>	4.2 mm x 11.5 mm Straight Drill
<b>15255K</b>	4.2 mm x 13 mm Straight Drill
<b>15101K</b>	1.8 mm/2.8 mm Step Drill
<b>15256K</b>	2.8 mm/3.3 mm Step Drill
<b>15257K</b>	3.3 mm/4.2 mm Step Drill
<b>15258K</b>	3.3 mm Countersink Drill
<b>15259K</b>	4.0 mm Countersink Drill
<b>15260K</b>	5.0 mm Countersink Drill
<b>15262K</b>	3.3 x 15 mm Straight Surgical Tap
<b>15264K</b>	4.0 x 15 mm Straight Surgical Tap
<b>15265K</b>	5.0 x 13 mm Straight Surgical Tap

#### Tissue Punch



#### Catalog #

#### Description

- 15237K** Ø3.5 mm, with Guide Pin
- 15238K** Ø4.1 mm, with Guide Pin
- 15239K** Ø5.0 mm, with Guide Pin

#### Round Marking Bur




#### Catalog #

#### Description

- L1405-23-00K** Ø2.3 mm

#### Initial Drills




#### Catalog #

#### Description

- 15200K** Ø1.5 x 15 mm, Non-Irrigated
- 15201K** Ø1.8 x 11.5 mm, Non-Irrigated
- 15202K** Ø1.8 x 15 mm, Non-Irrigated



**Parallel Pins**



*Catalog #*

*Description*

**15226K**  
**45045K**

Ø1.5/1.8 mm, Straight, 2-Pack  
Ø1.5/1.8 mm, 15°, 2-Pack

**Final Tapered Drills**



*Catalog #*

*Description*

**15203K**  
**15204K**

Ø3.0 x 13 mm  
Ø3.0 x 15 mm

**15205K**  
**15206K**  
**15207K**  
**15208K**

Ø3.5 x 10 mm  
Ø3.5 x 11.5 mm  
Ø3.5 x 13 mm  
Ø3.5 x 15 mm

**15209K**  
**15210K**  
**15211K**  
**15212K**

Ø4.1 x 10 mm  
Ø4.1 x 11.5 mm  
Ø4.1 x 13 mm  
Ø4.1 x 15 mm

**15213K**  
**15214K**  
**15215K**  
**15075K**

Ø5.0 x 10 mm  
Ø5.0 x 11.5 mm  
Ø5.0 x 13 mm  
Ø5.0 x 15 mm

**Removable Drill Stop  
for Tapered Drills**



*Catalog #*

*Description*

**15220K**  
**15221K**  
**15222K**  
**15223K**

3.0 mm  
3.5 mm  
4.1 mm  
5.0 mm

**Tapered Surgical Taps**



*Catalog #*

*Description*

**15217K**  
**15247K**

Ø3.5 x 10/11.5 mm  
Ø3.5 x 13/15 mm

**15218K**  
**15248K**

Ø4.1 x 10/11.5 mm  
Ø4.1 x 13/15 mm

**15219K**  
**15076K**

Ø5.0 x 10/11.5 mm  
Ø5.0 x 13/15 mm

Only PrimaConnex® Initial Drills are non-irrigated. All other Prima Tapered Drills, Tapered Taps, Straight Drills, Step Drills, Countersink Drills and Straight Taps are irrigated.





### Straight Drills



Catalog #	Description
15250K	Ø2.8 x 11.5 mm
15251K	Ø2.8 x 15 mm
15252K	Ø3.3 x 11.5 mm
15253K	Ø3.3 x 15 mm
15254K	Ø4.2 x 11.5 mm
15255K	Ø4.2 x 13 mm

### Step Drills for Straight Implants



Catalog #	Description
15101K	Ø1.8/Ø2.8 mm
15256K	Ø2.8/Ø3.3 mm
15257K	Ø3.3/Ø4.2 mm

### Countersink Drills for Straight Implants



Catalog #	Description
15258K	For Ø3.3 mm Straight Implants
15259K	For Ø4.0 mm Straight Implants
15260K	For Ø5.0 mm Straight Implants

### Straight Surgical Taps



Catalog #	Description
15261K	Ø3.3 x 11.5 mm
15262K	Ø3.3 x 15 mm
15263K	Ø4.0 x 11.5 mm
15264K	Ø4.0 x 15 mm
15265K	Ø5.0 x 13 mm

### Drill Extender



Catalog #
L1404-01K

### Tissue Punch



Catalog #	Description
15236K	Ø3.0 mm Tissue Punch with Guide Pin
15237K	Ø3.5 mm Tissue Punch with Guide Pin
15238K	Ø4.1 mm Tissue Punch with Guide Pin
15239K	Ø5.0 mm Tissue Punch with Guide Pin

### Bone Profiler



Catalog #	Description
15279K	Ø3.5 mm, with Guide Pin
15280K	Ø4.1 mm, with Guide Pin
15281K	Ø5.0 mm, with Guide Pin

Only PrimaConnex® Initial Drills are non-irrigated. All other Prima Tapered Drills, Tapered Taps, Straight Drills, Step Drills, Countersink Drills and Straight Taps are irrigated.



**PrimaConnex®  
Implant Driver, Ratchet**



*Catalog #*

*Description*

**15282K**  
**15283K**  
**15284K**

Ø3.5 mm  
Ø4.1 mm  
Ø5.0 mm

**PrimaConnex®  
Implant Driver, Latch**



*Catalog #*

*Description*


*Length*

**15286K**  
**15287K**  
**15288K**  
**15289K**  
**15290K**

Ø3.5/3.8 mm Implant Driver, Short  
Ø3.5/3.8 mm Implant Driver, Long  
Ø4.1/4.5 mm Implant Driver, Short  
Ø4.1/4.5 mm Implant Driver, Long  
Ø5.0/5.5/6.5/7.0/8.0/9.0 mm Implant Driver, Short

16 mm  
21 mm  
16 mm  
21 mm  
16 mm

**PrimaSolo®  
Implant Drivers**



*Catalog #*

*Description*

**15227K**  
**15228K**  
**15229K**  
**15230K**

Ø3.0 mm, Implant Driver  
Ø3.5 mm, Implant Driver  
Ø4.1 mm, Implant Driver  
Ø5.0 mm, Implant Driver

**Surgical Ratchet**



*Catalog #*

*Description*

**L1702-01K**

Surgical Ratchet

**Surgical Ratchet/  
Torque Wrench**



*Catalog #*

*Description*

**K15695**

Surgical Ratchet/Torque Wrench

**Quad Drivers**



*Catalog #*

*Description*

**15272K**  
**15273K**

Swivel Head, Short  
Swivel Head, Long

**15274K**  
**15275K**

Latch Type, Short  
Latch Type, Long

**15276K**  
**15277K**

Torque Wrench, Short  
Torque Wrench, Long



**Ratchet Adapter**

*Catalog #*

**15224K**



**Directional Indicators**

*Catalog #*

**45041K**

*Description*

Try-In Set




**Conversion Handles**

*Catalog #*

**L1407K**

*Description*

Conversion Handle, Latch Type



**Surgical Hand Driver**

*Catalog #*

**15271K**



**Depth Probe**

*Catalog #*


**15225K**



**Irrigation Port Cleaner**

*Catalog #*

**L1800-22-02K**



**Osteotomes**

*Catalog #*

**15063K**

**15067K**

**15068K**

**15069K**

**15070K**

*Description*

Osteotome Kit, Round, with Tray

1.8-2.5 mm, Round, Condensing

1.9-3.0 mm, Round, Condensing

2.2-3.6 mm, Round, Condensing

2.8-4.5 mm, Round, Condensing

**15064K**

**15071K**

**15072K**

**15073K**

**15074K**

Osteotome Kit, Sinus, with Tray

1.8-2.5 mm, Concave, Sinus

1.9-3.0 mm, Concave, Sinus

2.2-3.6 mm, Concave, Sinus


2.8-4.5 mm, Concave, Sinus

**15066K**

Osteotome Tray (empty)



## Surgical Kit Overview



<i>Catalog #</i>	<i>Description</i>
<b>15651K</b>	Tapered PrimaConnex® Guided Surgery Kit (with instruments) Includes the instrumentation required to place PrimaConnex® Tapered Implants
<b>15650K</b>	Guided Surgery Cassette (without instruments)

### Anchor Pins & Anchor Pin Drill



#### Catalog #

**E40100**  
**E20101**

#### Description

Anchor Pins (package of 3)  
Anchor Pin Drill Ø1.6 mm

### Tissue Punch



#### Catalog #

**15671K**  
**15672K**  
**15673K**

#### Description

PrimaConnex® Guided Tissue Punch, Ø3.5  
PrimaConnex® Guided Tissue Punch, Ø4.1  
PrimaConnex® Guided Tissue Punch, Ø5.0

### Guided Surgery Drills



#### Catalog #

**15674K**  
**15675K**

#### Description

PrimaConnex® Guided Initial Drill, Short  
PrimaConnex® Guided Initial Drill, Long

**15652K**  
**15653K**  
**15654K**  
**15655K**

PrimaConnex® Guided Surgery, Tapered Drill, Ø3.5 x 10  
PrimaConnex® Guided Surgery, Tapered Drill, Ø3.5 x 11.5  
PrimaConnex® Guided Surgery, Tapered Drill, Ø3.5 x 13  
PrimaConnex® Guided Surgery, Tapered Drill, Ø3.5 x 15

**15656K**  
**15657K**  
**15658K**  
**15659K**


PrimaConnex® Guided Surgery, Tapered Drill, Ø4.1 x 10  
PrimaConnex® Guided Surgery, Tapered Drill, Ø4.1 x 11.5  
PrimaConnex® Guided Surgery, Tapered Drill, Ø4.1 x 13  
PrimaConnex® Guided Surgery, Tapered Drill, Ø4.1 x 15

**15660K**  
**15661K**  
**15662K**  
**15663K**

PrimaConnex® Guided Surgery, Tapered Drill, Ø5.0 x 10  
PrimaConnex® Guided Surgery, Tapered Drill, Ø5.0 x 11.5  
PrimaConnex® Guided Surgery, Tapered Drill, Ø5.0 x 13  
PrimaConnex® Guided Surgery, Tapered Drill, Ø5.0 x 15



**Drill Guides**



*Catalog #*

*Description*


- 15676K**
- 15678K**
- 15679K**

- PrimaConnex® Guided Drill Guide, Initial Drill, Ø3.5
- PrimaConnex® Guided Drill Guide, Initial Drill, Ø4.1
- PrimaConnex® Guided Drill Guide, Initial Drill, Ø5.0

- 15680K**
- 15681K**
- 15682K**

- PrimaConnex® Guided Drill Guide, Ø3.5/Ø4.1
- PrimaConnex® Guided Drill Guide, Ø3.5/Ø5.0
- PrimaConnex® Guided Drill Guide, Ø4.1/Ø5.0

**Guided Surgery Bone Taps**



*Catalog #*

*Description*

- 15664K**
- 15665K**
- 15667K**
- 15668K**
- 15669K**
- 15670K**

- PrimaConnex® Guided Bone Tap, Ø3.5, Short
- PrimaConnex® Guided Bone Tap, Ø3.5, Long
- PrimaConnex® Guided Bone Tap, Ø4.1, Short
- PrimaConnex® Guided Bone Tap, Ø4.1, Long
- PrimaConnex® Guided Bone Tap, Ø5.0, Short
- PrimaConnex® Guided Bone Tap, Ø5.0, Long

**Guided Implant Drivers**



*Catalog #*

*Description*

- 15683K**
- 15684K**
- 15685K**

- PrimaConnex® Guided Implant Driver, Ø3.5
- PrimaConnex® Guided Implant Driver, Ø4.1
- PrimaConnex® Guided Implant Driver, Ø5.0

**Surgical Ratchet & Ratchet Adapter**




*Catalog #*

*Description*

- K15695**
- 15686K**

- Surgical Ratchet/Torque Wrench
- Ratchet Adapter

**Quad Drivers**



*Catalog #*

*Description*


- 15272K**
- 15273K**
- 15276K**
- 15277K**

- Quad Swivel, Short
- Quad Swivel, Long
- Quad Driver, Torque Wrench, Short
- Quad Driver, Torque Wrench, Long





Implantmed Surgical Motor & Accessories



<i>Catalog #</i>	<i>Description</i>
<b>WH-10207550</b>	WI-75 E/KM Contra-Angle 20:1 Handpiece for implant surgery (15-2,000 rpm)
<b>WH-30065000</b>	WS-75 Contra-Angle 20:1 Dismantlable Handpiece for implant surgery (15-2,000 rpm)
<b>WH-30032000</b>	WS-75 L G Contra-Angle 20:1 Dismantlable Handpiece with Mini LED+ for implant surgery (15-2,000 rpm)
<b>WH-30056000</b>	S-11 1:1 High Speed Straight Handpiece for oral surgery (300-40,000 rpm)
<b>WH-04363600</b>	Irrigation Tubing and Spike Assembly, 6-Pack, Disposable (for use with Implantmed series & Elcomed SA-310 motors)
<b>WH-00929300</b>	Irrigation Tubing with Luer Lock, Green (10-Pack)
<b>WH-04019000</b>	Irrigation Tubing Clips, 5-Pack (for use with Implantmed series motors to hold tubing to motor cable)
<b>WH-02610500</b>	Irrigation Needle/Canula (for use with internally irrigated drills)
<b>WH-06016400</b>	Y-Connector Complete
<b>WH-06946300</b>	Irrigation Spray Clip (new series), Left Side, E/KM, Green (3-Pack)
<b>WH-06948300</b>	Irrigation Spray Clip (new series), Right Side, E/KM, Green (3-Pack)
<b>WH-06948400</b>	Irrigation Spray Clip (new series), Left Side, Green (3-Pack)
<b>WH-06949500</b>	Irrigation Spray Clip (new series), Right Side, Green (3-Pack)
<b>WH-04757100</b>	Irrigation Spray Clip (old series), E/KM, Black (3-Pack)
<b>WH-05051000</b>	Irrigation Spray Clip (old series), Black (3-Pack)
<b>WH-02139800</b>	Bur Testing Gauge (to check all 2.35 mm & 1.6 mm burs)
<b>WH-10940011</b>	W&H F-1 Service Oil, 400 ml (Special shipping & handling required [aerosol can]. Contact Customer Service for shipping rates.)
<b>WH-02038200</b>	Spray Cap for Service Oil (reusable replacement)
<b>WH-04013500</b>	Sterilization Case

May not be available in all markets.





### Implantmed Surgical Motor & Accessories



<i>Catalog #</i>	<i>Description</i>
<b>WH-90000207</b>	SI-1015 Implantmed Base Motor Unit
<b>WH-90000213</b>	SI-1015 Implantmed Motor Unit Fully Loaded
<b>WH-30264001</b>	S-NW +SPI WiFi Foot Pedal
<b>WH-30058000</b>	S-11 L 1:1 Straight
<b>WH-30033000</b>	WS-75 L 20:1, Contra-angle

### Implantmed Surgical Motor & Accessories\*



<b>WH-90000208</b>	SI-1015 Implantmed Base Motor Unit w/Osstell
<b>WH-90000209</b>	SI-1015 Implantmed Motor Unit, WiFi Footpedal Upgrade
<b>WH-90000210</b>	SI-1015 Implantmed Motor Unit, WiFi Footpedal Upgrade & Osstell
<b>WH-90000211</b>	SI-1015 Implantmed Motor Unit, LED+ Motor Upgrade
<b>WH-90000212</b>	SI-1015 Implantmed Motor Unit, LED+ Motor Upgrade & Osstell
<b>WH-90000214</b>	SI-1015 Implantmed Motor Unit Fully Loaded & Osstell
<b>WH-30185000</b>	EM -19 Motor, 1.8m Cable, without LED+ Port
<b>WH-30281000</b>	EM-19 LC Motor, 1.8m Cable, LED+ Port
<b>WH-30264000</b>	S-NW + CAN WiFi Foot Pedal
<b>WH-30210000</b>	Osstell ISQ Module

May not be available in all markets.

\*Products available upon request.

## Terms and Conditions

**TERMS AND CONDITIONS OF SALE.** The following are the terms and conditions by which Keystone Dental, Inc. sells its products. As used herein, the term "Keystone" means Keystone Dental, Inc. and its applicable affiliates; the term "Products" means the Keystone products and services described in this invoice, and "Purchaser" means the person or company purchasing the Products.

**GOVERNING TERMS.** Any purchase and shipment of Products shall be deemed to be on the terms and conditions stated herein. Any and all terms and conditions submitted by Purchaser are hereby rejected.

**TAXES.** The prices set forth herein do not include any sales, use, excise, ad valorem, property or other taxes applicable to the sale, use or delivery of the Products, all of which shall be paid by Purchaser separately or added to the contract price and paid by Purchaser to Keystone.

**PRICES.** The prices set forth for the Products are in United States Dollars. All prices shall be as of the invoice date, regardless of the date the Product is originally ordered from Keystone.

**DELIVERY.** Products sold hereunder shall be delivered F.O.B. shipping point. Purchaser shall assume responsibility for all subsequent delivery or shipping charges.

**ACCEPTANCE.** Purchaser must inspect all Products for defects and/or non-conformity within seven (7) calendar days of delivery. After this period, Purchaser will have accepted such Products.

**PAYMENT.** Payment terms shall be as set forth on the front of the invoice. If payment is late, the maximum statutory interest rate will apply on such late amount. Keystone may, in its sole discretion, suspend shipments of Products until full payment by Purchaser of all past due amounts to Keystone is made.

**RETURN; EXCHANGE.** Purchaser may return any Product within thirty (30) calendar days of the invoice date so long as there is an offsetting (exchange) order of equal or greater value. Purchaser may return any Product after thirty (30) calendar days from the invoice date, but no later than ninety (90) calendar days after the invoice date, so long as there is an offsetting (exchange) order of equal or greater value, provided, however, such returns shall be subject to a 15% re-stocking fee. Keystone will not accept any returns or exchanges after ninety (90) calendar days following the invoice date.

No biomaterial products (bone grafting and membrane) will be accepted for return or exchange. No returns or exchanges will be accepted for Products that (i) were originally sold as part of a kit, pack or set, (ii) were included on the original invoice at no charge, (iii) have less than six (6) months remaining before expiration or (iv) were custom ordered. The Purchaser shall pay all shipping costs associated with the return or exchange of a Product. No Products will be accepted by Keystone for return or exchange without an RMA number given to Purchaser by Keystone Customer Service. Returned Products must be unexpired, unopened and undamaged, returned freight prepaid.

**WARRANTY.** EXCEPT AS EXPRESSLY SET FORTH IN KEYSTONE'S LIFETIME IMPLANT WARRANTY PROGRAM, KEYSTONE DOES NOT MAKE ANY WARRANTY, EXPRESS OR IMPLIED, EXCEPT THAT ITS PRODUCTS SHALL BE FREE FROM DEFECTS IN MATERIAL AND/OR WORKMANSHIP. THIS WARRANTY APPLIES ONLY TO THE ORIGINAL PURCHASER. PURCHASER ASSUMES ALL RISKS AND LIABILITY RESULTING FROM THE USE OF THESE PRODUCTS, WHETHER USED SEPARATELY OR IN COMBINATION WITH OTHER PRODUCTS. THE WARRANTY DOES NOT APPLY IF THE FAILURE OF THE PRODUCT IS CAUSED BY TRAUMA OR ACCIDENT OR BY THE PATIENT IN WHOM THE IMPLANT WAS INSTALLED OR IF THE PURCHASER IS NOT CURRENT ON ALL OUTSTANDING PAYMENTS TO KEYSTONE AT THE TIME OF SUCH CLAIM. KEYSTONE MAY MODIFY OR TERMINATE ITS WARRANTY PROGRAM AT ANY TIME, HOWEVER ANY CHANGES TO THE WARRANTY PROGRAM WILL NOT AFFECT PRODUCTS ALREADY PURCHASED AT THE TIME OF SUCH CHANGE. PURCHASER'S SOLE REMEDY WITH RESPECT TO CLAIMS RELATED TO KEYSTONE PRODUCTS IS RESTRICTED TO THE WARRANTY POLICY.

**LIMITATION OF LIABILITY.** UNDER NO CIRCUMSTANCES SHALL KEYSTONE OR ANY OF ITS AFFILIATES BE LIABLE FOR BUSINESS INTERRUPTION, LOSS OF PROFITS, SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE AND FROM ANY CAUSE WHETHER BASED IN CONTRACT, WARRANTY OR TORT (INCLUDING NEGLIGENCE) OR OTHER LEGAL THEORY, EVEN IF KEYSTONE HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

**INTELLECTUAL PROPERTY.** All intellectual property rights in Products sold to Purchaser shall remain the exclusive property of Keystone and its licensors.

**PERSONAL DATA.** It is Purchaser's responsibility to obtain patient's consent to transfer any personal data to Keystone.

**TERMINATION.** Keystone may terminate this or any agreement with Purchaser to purchase its products upon immediate written notice if Purchaser fails to pay any of its debt obligations to Keystone on time.

**FORCE MAJEURE.** Keystone shall not be liable for any loss or damage due to failure or delay arising out of any cause beyond the reasonable control of Keystone. In the event of any failure or delay resulting from such causes, an equitable adjustment of delivery and any other appropriate terms and conditions shall be made. No such failure or delay shall be the basis for an increase in Keystone's obligations, nor any termination by Purchaser.

**VALIDITY.** If any provision of these Terms and Conditions is found to be illegal or unenforceable in any respect, such illegality or unenforceability shall not affect any other provision of these Terms and Conditions, all of which shall remain enforceable in accordance with their terms.

**GENERAL PROVISIONS.** The purchase of the Products and these Terms and Conditions are governed by the laws of the Commonwealth of Massachusetts, without regard to conflicts of law principles. Any dispute, claim or controversy arising out of or in connection with the purchase of the Products shall be resolved by arbitration in Boston, Massachusetts. This document incorporates all oral and written representations between the parties and constitutes the entire agreement and understanding of the parties with respect to the subject matter hereof and supersedes any and all other agreements either oral or written between the parties with respect to such subject matter. No amendment or modification may be made to this document unless in writing and duly executed by an authorized representative of Keystone. Keystone may subcontract, assign or transfer its obligations or rights to a competent third party. Purchaser may not assign or transfer any of its obligations without the consent of Keystone.



**Keystone Dental, Inc.**  
**Global Headquarters**

154 Middlesex Turnpike  
Burlington, MA 01803 USA  
General Inquiries: 781-328-3600  
Customer Relations: 866-902-9272 or 781-328-3490  
Fax: 866-903-9272 or 781-328-3400  
Email: info@keystonedental.com



Caution, consult accompanying documents

**Trademark Acknowledgements**

TiLobe, PrimaPlus, Quick-Cap and "Smarter Thinking. Simpler Design." are trademarks and registered trademarks of Keystone Dental, Inc.